

P62/SQSTM1 Monoclonal antibody

【货号】 AM20109

【规格】 50ul

【产品信息】

◆产品名称	P62/SQSTM1 Monoclonal antibody
◆规格	50ul
◆别名	A170, EBIAP, ORCA, OSIL, p60, P62, P62/SQSTM1, P62; SQSTM1, p62B, PDB3, sequestosome 1, SQSTM1, Ubiquitin binding protein p62, ZIP3
◆产品类型	Monoclonal
◆基因全称	sequestosome 1
◆经测试应用	WB, IF, IHC
◆宿主	Mouse
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human P62/SQSTM1
◆理论分子量	48 kDa
◆实际分子量	55-60 kDa
◆Gene ID (NCBI)	8878
◆GenBank Accession Number	BC017222
◆UNIPROT	Q13501
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:5000-1:50000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF) 1:2000-1:8000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。