

E-cadherin Monoclonal antibody

【货号】 AM20111

【规格】 50ul

【产品信息】

◆产品名称	E-cadherin Monoclonal antibody
◆规格	50ul
◆别名	Arc 1; Cadherin 1; CD324; CDH1; CDHE; E cadherin; ECAD; E-cad; E-cadherin
◆产品类型	Monoclonal
◆基因全称	cadherin 1, type 1, E-cadherin (epithelial)
◆经测试应用	WB, IF, IHC
◆宿主	Mouse
◆种属反应性	Human, Rat
◆免疫原	A synthesized peptide derived from human E-cadherin
◆理论分子量	97 kDa
◆实际分子量	120 kDa
◆Gene ID (NCBI)	999
◆GenBank Accession Number	BC141838
◆UNIPROT	P12830
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2000-1:8000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF) 1:1000-1:4000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。