

HO-1/HMOX1 Monoclonal antibody

【货号】 AM20127

【规格】 50ul/100ul

【产品信息】

◆产品名称	HO-1/HMOX1 Monoclonal antibody
◆规格	50ul/100ul
◆别名	HMOX1, HO-1, Heme oxygenase 1 soluble form, Heme oxygenase 1, heme oxygenase (decycling) 1
◆产品类型	Monoclonal
◆基因全称	heme oxygenase (decycling) 1
◆经测试应用	WB, IHC, IF-P
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat, Pig
◆免疫原	A synthesized peptide derived from human HO-1
◆理论分子量	33 kDa
◆实际分子量	33 kDa
◆Gene ID (NCBI)	3162
◆GenBank Accession Number	BC001491
◆UNIPROT	P09601
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:250-1:1000 for Immunohistochemistry (IHC) 1:100-1:400 for Immunofluorescence (IF-P)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。