

Arginase-1 Monoclonal antibody

【货号】 AM20135

【规格】 50ul/100ul

【产品信息】

◆产品名称	Arginase-1 Monoclonal antibody
◆规格	50ul/100ul
◆别名	ARG1, Arginase 1, liver Arginase, arginase, liver, Arg-1
◆产品类型	Monoclonal
◆基因全称	arginase, liver
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human Arginase-1
◆理论分子量	35 kDa
◆实际分子量	35 kDa
◆Gene ID (NCBI)	383
◆GenBank Accession Number	BC005321
◆UNIPROT	P05089
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:100-1:400 for Immunofluorescence (IF/P) 1:1000-1:2500 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。