

GLUT1 Monoclonal antibody

【货号】 AM20151

【规格】 50ul/100ul

【产品信息】

◆产品名称	GLUT1 Monoclonal antibody
◆规格	50ul/100ul
◆别名	SLC2A1, SLC2A1, GLUT1, GLUT-1, 2A5A2
◆产品类型	Monoclonal
◆基因全称	solute carrier family 2 (facilitated glucose transporter), member 1
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Mouse
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human GLUT1
◆理论分子量	54 kDa
◆实际分子量	50-70 kDa
◆Gene ID (NCBI)	6513
◆GenBank Accession Number	BC121804
◆UNIPROT	P11166
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:500-1:1500 for western blotting (WB) 1:500-1:2000 for Immunohistochemistry (IHC) 1:100-1:400 for Immunofluorescence (IF)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。