

HSP70 Monoclonal antibody

【货号】 AM20165

【规格】 50ul/100ul

【产品信息】

◆产品名称	HSP70 Monoclonal antibody
◆规格	50ul/100ul
◆别名	HSPA1A, 4E7E5, Heat shock 70 kDa protein 1A, heat shock 70kDa protein 1A, Heat shock
◆产品类型	Monoclonal
◆基因全称	heat shock 70kDa protein 1A
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human HSP70
◆理论分子量	70 kDa
◆实际分子量	60-70 kDa
◆Gene ID (NCBI)	3303
◆GenBank Accession Number	BC009322
◆UNIPROT	P0DMV8
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:10000-1:50000 for western blotting (WB) 1:200-1:800 for Immunohistochemistry (IHC) 1:200-1:800 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。