

## AGER/RAGE Monoclonal antibody

【货号】 AM20172

【规格】 50ul/100ul

### 【产品信息】

◆产品名称	AGER/RAGE Monoclonal antibody
◆规格	50ul/100ul
◆别名	AGER, RAGE, 3D4G11, Advanced glycosylation end product-specific receptor
◆产品类型	Monoclonal
◆基因全称	advanced glycosylation end product-specific receptor
◆经测试应用	WB, IHC, IF-P
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human AGER
◆理论分子量	43 kDa
◆实际分子量	45-55 kDa
◆Gene ID (NCBI)	177
◆GenBank Accession Number	BC020669
◆UNIPROT	Q15109
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:500-1:3000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF-P) 1:75-1:300 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。