

MLKL Monoclonal antibody

【货号】 AM20179

【规格】 50ul/100ul

【产品信息】

◆产品名称	MLKL Monoclonal antibody
◆规格	50ul/100ul
◆别名	3D4C6, hMLKL, Mixed lineage kinase domain-like protein
◆产品类型	Monoclonal
◆基因全称	mixed lineage kinase domain-like
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Mouse
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human MLKL
◆理论分子量	54 kDa
◆实际分子量	45-55 kDa
◆Gene ID (NCBI)	197259
◆GenBank Accession Number	BC028141
◆UNIPROT	Q8NB16
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2000-1:20000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunofluorescence (IF-P) 1:250-1:1000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。