

YAP1 Monoclonal antibody

【货号】 AM20193

【规格】 50ul/100ul

【产品信息】

◆产品名称	YAP1 Monoclonal antibody
◆规格	50ul/100ul
◆别名	yes associated protein 1, YAP, Transcriptional coactivator YAP1, Protein yorkie homolog, 65 kDa Yes associated protein
◆产品类型	Monoclonal
◆基因全称	Yes-associated protein 1, 65kDa
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human YAP1
◆理论分子量	54 kDa
◆实际分子量	70-75 kDa
◆Gene ID (NCBI)	10413
◆GenBank Accession Number	BC038235
◆UNIPROT	P46937
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:10000 for western blotting (WB) 1:500-1:2000 for Immunohistochemistry (IHC) 1:200-1:800 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。