

NLRP3 Monoclonal antibody

【货号】 AM20203

【规格】 50ul/100ul

【产品信息】

◆产品名称	NLRP3 Monoclonal antibody
◆规格	50ul/100ul
◆别名	3H1A7, Angiotensin/vasopressin receptor AII/AVP-like, C1orf7, CIAS1, Cold-induced autoinflammatory syndrome 1 protein
◆产品类型	Monoclonal
◆基因全称	NLR family, pyrin domain containing 3
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human NLRP3
◆理论分子量	118 kDa
◆实际分子量	110 kDa
◆Gene ID (NCBI)	114548
◆GenBank Accession Number	NM_001079821
◆UNIPROT	Q96P20
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:100-1:400 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。