

ICAM-1/CD54 Mouse Monoclonal antibody

【货号】 AM20263

【规格】 50ul/100ul

【产品信息】

◆产品名称	ICAM-1/CD54 Mouse Monoclonal antibody
◆规格	50ul/100ul
◆别名	CD54, ICAM-1, 2F9A8, BB2, ICAM 1
◆产品类型	Monoclonal
◆基因全称	intercellular adhesion molecule 1
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Mouse
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human ICAM-1
◆理论分子量	90 kDa
◆实际分子量	90-100 kDa
◆Gene ID (NCBI)	3383
◆GenBank Accession Number	BC015969
◆UNIPROT	P05362
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:100-1:400 for Immunofluorescence (IF-P)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。