

PI3 Kinase p85 Alpha Mouse Monoclonal antibody

【货号】 AM20267

【规格】 50ul/100ul

【产品信息】

◆产品名称	PI3 Kinase p85 Alpha Mouse Monoclonal antibody
◆规格	50ul/100ul
◆别名	P85a; PIK3R1, PI3K p85 alpha, PI3K p85 α , PIK3R1, 4G3C11
◆产品类型	Monoclonal
◆基因全称	phosphoinositide-3-kinase, regulatory subunit 1 (alpha)
◆经测试应用	WB, IHC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human PI3 Kinase p85 Alpha
◆理论分子量	85 kDa
◆实际分子量	85 kDa
◆Gene ID (NCBI)	5295
◆GenBank Accession Number	BC030815
◆UNIPROT	P27986
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:500-1:2000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。