

# LC3 Polyclonal antibody

【货号】 AP20104

【规格】 50ul

## 【产品信息】

◆产品名称	LC3 Polyclonal antibody
◆规格	50ul
◆别名	MAP1LC3B, MAP1A/MAP1B LC3 B, MAP1A/1BLC3, MAP1 light chain 3-like protein 2, LC3B
◆产品类型	Polyclonal
◆基因全称	microtubule-associated protein 1 light chain 3 beta
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human LC3
◆理论分子量	15 kDa
◆实际分子量	15 kDa
◆Gene ID (NCBI)	81631
◆GenBank Accession Number	BC067797
◆UNIPROT	Q9GZQ8
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:250-1:1000 for Immunofluorescence (IF) 1:50-1:500 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。