

## ZO-1 Polyclonal antibody

【货号】 AP20106

【规格】 50ul

### 【产品信息】

◆产品名称	ZO-1 Polyclonal antibody
◆规格	50ul
◆别名	TJP1, ZO1, ZO 1, Tight junction protein ZO-1, Tight junction protein ZO 1
◆产品类型	Polyclonal
◆基因全称	tight junction protein 1 (zona occludens 1)
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human ZO-1
◆理论分子量	195 kDa
◆实际分子量	230 kDa
◆Gene ID (NCBI)	7082
◆GenBank Accession Number	BC111712
◆UNIPROT	Q07157
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:5000-1:50000 for western blotting (WB) 1:100-1:2000 for Immunofluorescence (IF) 1:1000-1:4000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。