

Histone H3 Polyclonal antibody

【货号】 AP20118

【规格】 50ul/100ul

【产品信息】

◆产品名称	Histone H3 Polyclonal antibody
◆规格	50ul/100ul
◆别名	HIST2H3A, H3-clustered histone 15, H3-clustered histone 14, H3-clustered histone 13, H3C15
◆产品类型	Polyclonal
◆基因全称	histone cluster 2, H3a
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human Histone H3
◆理论分子量	15 kDa
◆实际分子量	15 kDa
◆Gene ID (NCBI)	333932
◆GenBank Accession Number	BC015544
◆UNIPROT	Q71DI3
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2000-1:16000 for western blotting (WB) 1:600-1:2400 for Immunofluorescence (IF/ICC) 1:50-1:500 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。