

# P16-INK4A Polyclonal antibody

【货号】 AP20120

【规格】 50ul/100ul

## 【产品信息】

|                              |  |
|------------------------------|--|
| ◆产品名称                        | P16-INK4A Polyclonal antibody  |
| ◆规格                          | 50ul/100ul   |
| ◆别名                          | p16, p16 INK4, P16 INK4A, p16-INK4, p16INK4A   |
| ◆产品类型                        | Polyclonal   |
| ◆基因全称                        | cyclin-dependent kinase inhibitor 2A   |
| ◆经测试应用                       | WB, IHC, FC  |
| ◆宿主                          | Rabbit   |
| ◆种属反应性                       | Human  |
| ◆免疫原                         | A synthesized peptide derived from human P16-INK4A   |
| ◆理论分子量                       | 16 kDa   |
| ◆实际分子量                       | 16 kDa   |
| ◆Gene ID (NCBI)              | 1029   |
| ◆GenBank<br>Accession Number | BC021998   |
| ◆UNIPROT                     | P42771   |
| ◆形式                          | Liquid   |
| ◆保存条件                        | Store at -20°C. Stable for one year after shipment.  |
| ◆建议稀释度                       | 1:1000-1:3000 for western blotting (WB)<br>1:500-1:2000 for Immunohistochemistry (IHC)<br>0.40 ug per 10 <sup>6</sup> cells in a 100 μl suspension for Flow Cytometry (FC) (INTRA) |

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。