

## SLC7A11/xCT Polyclonal antibody

【货号】 AP20121

【规格】 50ul/100ul

### 【产品信息】

◆产品名称	SLC7A11/xCT Polyclonal antibody
◆规格	50ul/100ul
◆别名	SLC7A11, xCT, Amino acid transport system xc, Amino acid transport system xc-, CCBR1
◆产品类型	Polyclonal
◆基因全称	solute carrier family 7, (cationic amino acid transporter, y+ system) member 11
◆经测试应用	WB, IHC, IP
◆宿主	Rabbit
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human SLC7A11
◆理论分子量	55 kDa
◆实际分子量	45-55 kDa
◆Gene ID (NCBI)	23657
◆GenBank Accession Number	BC012087
◆UNIPROT	Q9UPY5
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:500-1:2000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。