

iNOS Polyclonal antibody

【货号】 AP20130

【规格】 50ul/100ul

【产品信息】

◆产品名称	iNOS Polyclonal antibody
◆规格	50ul/100ul
◆别名	NOS2, EC:1.14.13.39, Hepatocyte NOS, HEP-NOS, Inducible NO synthase
◆产品类型	Polyclonal
◆基因全称	nitric oxide synthase 2, inducible
◆经测试应用	WB, IF/ICC,FC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human iNOS
◆理论分子量	131 kDa
◆实际分子量	110-130 kDa , 60-70 kDa
◆Gene ID (NCBI)	4843
◆GenBank Accession Number	NM_000625
◆UNIPROT	P35228
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:250-1:1000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF/ICC) 0.40 ug per 10 ⁶ cells in a 100 μl suspension for Flow Cytometry (FC) (INTRA)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。