

LDHA-Specific Polyclonal antibody

【货号】 AP20131

【规格】 50ul/100ul

【产品信息】

◆产品名称	LDHA-Specific Polyclonal antibody
◆规格	50ul/100ul
◆别名	LDHA, Cell proliferation-inducing gene 19 protein, EC:1.1.1.27, lactate dehydrogenase A, LDH A
◆产品类型	Polyclonal
◆基因全称	lactate dehydrogenase A
◆经测试应用	WB, IF/ICC,IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human LDHA
◆理论分子量	37 kDa
◆实际分子量	30-37 kDa
◆Gene ID (NCBI)	3939
◆GenBank Accession Number	NM_005566
◆UNIPROT	P00338
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:10-1:100 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。