

Beclin 1 Polyclonal antibody

【货号】 AP20137

【规格】 50ul/100ul

【产品信息】

◆产品名称	Beclin 1 Polyclonal antibody
◆规格	50ul/100ul
◆别名	BECN1, Beclin-1-C 37 kDa, Beclin-1-C 35 kDa, Beclin-1, Beclin1
◆产品类型	Polyclonal
◆基因全称	beclin 1, autophagy related
◆经测试应用	WB, IF/ICC, IF-P, IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human Beclin 1
◆理论分子量	52 kDa
◆实际分子量	56-60 kDa
◆Gene ID (NCBI)	8678
◆GenBank Accession Number	BC010276
◆UNIPROT	Q14457
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:500-1:5000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF-P) 1:100-1:400 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。