

LIAS Polyclonal antibody

【货号】 AP20138

【规格】 50ul/100ul

【产品信息】

| | |
|------------------------------|--|
| ◆产品名称 | LIAS Polyclonal antibody |
| ◆规格 | 50ul/100ul |
| ◆别名 | HUSSY 01, LAS, LIAS, Lip syn, LIP1, Lipoate synthase, Lipoic acid synthase, lipoic acid synthetase, Lipoyl synthase, mitochondrial, LS |
| ◆产品类型 | Polyclonal |
| ◆基因全称 | lipoic acid synthetase |
| ◆经测试应用 | WB, IHC |
| ◆宿主 | Rabbit |
| ◆种属反应性 | Human, Mouse, Rat |
| ◆免疫原 | A synthesized peptide derived from human LIAS |
| ◆理论分子量 | 42 kDa |
| ◆实际分子量 | 34-42 kDa |
| ◆Gene ID (NCBI) | 11019 |
| ◆GenBank Accession Number | BC023635 |
| ◆UNIPROT | O43766 |
| ◆形式 | Liquid |
| ◆保存条件 | Store at -20℃. Stable for one year after shipment. |
| ◆建议稀释度 | 1:500-1:2500 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) |

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。