

# MYH10-Specific Polyclonal antibody

【货号】 AP20161

【规格】 50ul/100ul

## 【产品信息】

◆产品名称	MYH10-Specific Polyclonal antibody
◆规格	50ul/100ul
◆别名	Myosin 10, MYH10, NMMHC II-b, Myosin-10, Myosin heavy chain, non-muscle IIb
◆产品类型	Polyclonal
◆基因全称	myosin, heavy chain 10, non-muscle
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human MYH10
◆理论分子量	231 kDa
◆实际分子量	231 kDa
◆Gene ID (NCBI)	4628
◆GenBank Accession Number	NM_005964
◆UNIPROT	P35580
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:5000-1:50000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:10-1:50 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。