

XBP1S-specific Polyclonal antibody

【货号】 AP20164

【规格】 50ul/100ul

【产品信息】

◆产品名称	XBP1S-specific Polyclonal antibody
◆规格	50ul/100ul
◆别名	XBP1, Tax-responsive element-binding protein 5, TREB-5, X box binding protein 1, X-box-binding protein 1
◆产品类型	Polyclonal
◆基因全称	X-box binding protein 1
◆经测试应用	WB, IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human XBP1
◆理论分子量	29 kDa
◆实际分子量	55 kDa
◆Gene ID (NCBI)	7494
◆GenBank Accession Number	BC000938
◆UNIPROT	P17861
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:500-1:3000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。