

IGF2BP1 Polyclonal antibody

【货号】 AP20170

【规格】 50ul/100ul

【产品信息】

◆产品名称	IGF2BP1 Polyclonal antibody
◆规格	50ul/100ul
◆别名	IMP 1, IGF-II mRNA-binding protein 1, IGF2BP 1, IGF2 mRNA-binding protein 1, CRD-BP
◆产品类型	Polyclonal
◆基因全称	insulin-like growth factor 2 mRNA binding protein 1
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human IGF2BP1
◆理论分子量	65 kDa
◆实际分子量	63-70 kDa
◆Gene ID (NCBI)	10642
◆GenBank Accession Number	BC156957
◆UNIPROT	Q9NZI8
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:100-1:400 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。