

# MCT1 Polyclonal antibody

【货号】 AP20196

【规格】 50ul/100ul

## 【产品信息】

◆产品名称	MCT1 Polyclonal antibody
◆规格	50ul/100ul
◆别名	SLC16A1, FLJ36745, HHF7, MCT, MCT 1
◆产品类型	Polyclonal
◆基因全称	solute carrier family 16, member 1 (monocarboxylic acid transporter 1)
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human MCT1
◆理论分子量	54 kDa
◆实际分子量	34-42 kDa
◆Gene ID (NCBI)	6566
◆GenBank Accession Number	BC026317
◆UNIPROT	P53985
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:6000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF/ICC) 1:400-1:1600 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。