

# MMP13 Polyclonal antibody

【货号】 AP20197

【规格】 50ul/100ul

## 【产品信息】

◆产品名称	MMP13 Polyclonal antibody
◆规格	50ul/100ul
◆别名	CLG3, Collagenase 3, EC:3.4.24.-, Matrix metalloproteinase 13, Matrix metalloproteinase-13
◆产品类型	Polyclonal
◆基因全称	matrix metalloproteinase 13
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human MMP13
◆理论分子量	54 kDa
◆实际分子量	65-70 kDa, 45-55 kDa, 40 kDa
◆Gene ID (NCBI)	4322
◆GenBank Accession Number	BC074808
◆UNIPROT	P45452
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:500-1:2000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。