

SCD1 Polyclonal antibody

【货号】 AP20202

【规格】 50ul/100ul

【产品信息】

◆产品名称	SCD1 Polyclonal antibody
◆规格	50ul/100ul
◆别名	SCD, EC:1.14.19.1, Delta-9 desaturase, Delta(9)-desaturase, Delta(9) desaturase
◆产品类型	Polyclonal
◆基因全称	stearoyl-CoA desaturase (delta-9-desaturase)
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human SCD
◆理论分子量	41 kDa
◆实际分子量	35-42 kDa
◆Gene ID (NCBI)	6319
◆GenBank Accession Number	BC005807
◆UNIPROT	O00767
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:8000 for western blotting (WB) 1:100-1:400 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。