

GBP4 Polyclonal antibody

【货号】 AP20205

【规格】 50ul/100ul

【产品信息】

◆产品名称	GBP4 Polyclonal antibody
◆规格	50ul/100ul
◆别名	EC:3.6.5.-, GBP-4, GTP-binding protein 4, Guanine nucleotide-binding protein 4, Guanylate-binding protein 4
◆产品类型	Polyclonal
◆基因全称	guanylate binding protein 4
◆经测试应用	WB
◆宿主	Rabbit
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human GBP4
◆理论分子量	73 kDa
◆实际分子量	73 kDa
◆Gene ID (NCBI)	115361
◆GenBank Accession Number	BC070055
◆UNIPROT	Q96PP9
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:500-1:1500 for western blotting (WB)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。