

NR4A1 Polyclonal antibody

【货号】 AP20209

【规格】 50ul/100ul

【产品信息】

| | |
|------------------------------|---|
| ◆产品名称 | NR4A1 Polyclonal antibody |
| ◆规格 | 50ul/100ul |
| ◆别名 | Early response protein NAK1, GFRP1, HMR, N10, NAK 1, NAK1, NGFIB, NR4A1, NUR77, Orphan nuclear receptor HMR, Orphan nuclear receptor TR3, ST 59, Testicular receptor 3, TR3 |
| ◆产品类型 | Polyclonal |
| ◆基因全称 | nuclear receptor subfamily 4, group A, member 1 |
| ◆经测试应用 | WB, IHC |
| ◆宿主 | Rabbit |
| ◆种属反应性 | Human, Mouse |
| ◆免疫原 | A synthesized peptide derived from human NR4A1 |
| ◆理论分子量 | 64 kDa |
| ◆实际分子量 | 64-72 kDa |
| ◆Gene ID (NCBI) | 3164 |
| ◆GenBank Accession Number | BC016147 |
| ◆UNIPROT | P22736 |
| ◆形式 | Liquid |
| ◆保存条件 | Store at -20℃. Stable for one year after shipment. |
| ◆建议稀释度 | 1:250-1:500 for western blotting (WB) |

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。