

PPP1CC Polyclonal antibody

【货号】 AP20211

【规格】 50ul/100ul

【产品信息】

◆产品名称	PPP1CC Polyclonal antibody
◆规格	50ul/100ul
◆别名	EC:3.1.3.16, PP-1G, Serine/threonine-protein phosphatase PP1-gamma catalytic subunit
◆产品类型	Polyclonal
◆基因全称	protein phosphatase 1, catalytic subunit, gamma isoform
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human PPP1CC
◆理论分子量	35 kDa
◆实际分子量	35 kDa
◆Gene ID (NCBI)	5501
◆GenBank Accession Number	BC014073
◆UNIPROT	P36873
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:1000-1:8000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。