

# HLA-E Polyclonal antibody

【货号】 AP20222

【规格】 50ul/100ul

## 【产品信息】

◆产品名称	HLA-E Polyclonal antibody
◆规格	50ul/100ul
◆别名	HLA 6.2, HLA E, HLA-6.2, HLAE, MHC class I antigen E
◆产品类型	Polyclonal
◆基因全称	major histocompatibility complex, class I, E
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human HLA-E
◆理论分子量	40 kDa
◆实际分子量	40 kDa
◆Gene ID (NCBI)	3133
◆GenBank Accession Number	BC002578
◆UNIPROT	P13747
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:8000 for western blotting (WB) 1:250-1:1000 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。