

JNK Polyclonal antibody

【货号】 AP20225

【规格】 50ul/100ul

【产品信息】

◆产品名称	JNK Polyclonal antibody
◆规格	50ul/100ul
◆别名	c Jun N terminal kinase 1, JNK 46, JNK1, JNK1A2, JNK21B1/2, MAP kinase 8, MAPK 8, MAPK8, PRKM8, SAPK1
◆产品类型	Polyclonal
◆基因全称	mitogen-activated protein kinase 8
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human JNK
◆理论分子量	48 kDa
◆实际分子量	40-55 kDa
◆Gene ID (NCBI)	5599
◆GenBank Accession Number	NM_138982
◆UNIPROT	P45983
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:250-1:1500 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。