

# DHX9 Polyclonal antibody

【货号】 AP20240

【规格】 50ul/100ul

## 【产品信息】

◆产品名称	DHX9 Polyclonal antibody
◆规格	50ul/100ul
◆别名	ATP-dependent RNA helicase A, DDX9, DExH-box helicase 9, EC:3.6.4.13, LKP
◆产品类型	Polyclonal
◆基因全称	DEAH (Asp-Glu-Ala-His) box polypeptide 9
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human DHX9
◆理论分子量	141 kDa
◆实际分子量	141 kDa
◆Gene ID (NCBI)	1660
◆GenBank Accession Number	BC014246
◆UNIPROT	Q08211
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:125-1:500 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。