

METTL14 Recombinant antibody

【货号】 RM20141

【规格】 50ul/100ul

【产品信息】

◆产品名称	METTL14 Recombinant antibody
◆规格	50ul/100ul
◆别名	2E7, hMETTL14, Methyltransferase-like protein 14, METT L14, N6-adenosine-methyltransferase non-catalytic subunit
◆产品类型	Monoclonal
◆基因全称	methyltransferase like 14
◆经测试应用	WB, IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human METTL14
◆理论分子量	52 kDa
◆实际分子量	60-70 kDa
◆Gene ID (NCBI)	57721
◆GenBank Accession Number	BC007449
◆UNIPROT	Q9HCE5
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:7000 for western blotting (WB) 1:500-1:2000 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；

◆本产品仅限科研使用。