

ACSL4/FACL4 Recombinant antibody

【货号】 RM20147

【规格】 50ul/100ul

【产品信息】

◆产品名称	ACSL4/FACL4 Recombinant antibody
◆规格	50ul/100ul
◆别名	ACSL4, FACL4, EC:6.2.1.15, Arachidonate--CoA ligase, ACSL 4
◆产品类型	Monoclonal
◆基因全称	acyl-CoA synthetase long-chain family member 4
◆经测试应用	WB, IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human ACSL4
◆理论分子量	79 kDa
◆实际分子量	72 kDa
◆Gene ID (NCBI)	2182
◆GenBank Accession Number	BC034959
◆UNIPROT	O60488
◆形式	Liquid
◆保存条件	Store at -20°C . Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。