

IGF2BP2 Recombinant antibody

【货号】 RM20169

【规格】 50ul/100ul

【产品信息】

| | |
|------------------------------|--|
| ◆产品名称 | IGF2BP2 Recombinant antibody |
| ◆规格 | 50ul/100ul |
| ◆别名 | insulin like growth factor 2 mRNA binding protein 2, IMP-2, IMP2, IGF-II mRNA-binding protein 2, IGF2 mRNA-binding protein 2 |
| ◆产品类型 | Monoclonal |
| ◆基因全称 | insulin-like growth factor 2 mRNA binding protein 2 |
| ◆经测试应用 | WB, IHC |
| ◆宿主 | Rabbit |
| ◆种属反应性 | Human, Mouse, Rat |
| ◆免疫原 | A synthesized peptide derived from human IGF2BP2 |
| ◆理论分子量 | 65 kDa |
| ◆实际分子量 | 60-70 kDa |
| ◆Gene ID (NCBI) | 10644 |
| ◆GenBank Accession Number | BC021290 |
| ◆UNIPROT | Q9Y6M1 |
| ◆形式 | Liquid |
| ◆保存条件 | Store at -20°C. Stable for one year after shipment. |
| ◆建议稀释度 | 1:1000-1:5000 for western blotting (WB) 1:250-1:1000 for Immunohistochemistry (IHC) |

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。