

## HMGB1 Recombinant antibody

【货号】 RM20189

【规格】 50ul/100ul

### 【产品信息】

◆产品名称	HMGB1 Recombinant antibody
◆规格	50ul/100ul
◆别名	HMGB1, HMGB 1, High mobility group protein B1, High mobility group protein 1, high mobility group box 1
◆产品类型	Monoclonal
◆基因全称	high-mobility group box 1
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human HMGB1
◆理论分子量	25 kDa
◆实际分子量	25-30 kDa
◆Gene ID (NCBI)	3146
◆GenBank Accession Number	BC003378
◆UNIPROT	P09429
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:50-1:250 for Immunofluorescence (IF/ICC) 1:500-1:2000 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。