

ATF4 Recombinant antibody

【货号】 RM20199

【规格】 50ul/100ul

【产品信息】

◆产品名称	ATF4 Recombinant antibody
◆规格	50ul/100ul
◆别名	5N3, cAMP-dependent transcription factor ATF-4, CREB2, Cyclic AMP-responsive element-binding protein 2, Tax-responsive enhancer element-binding protein 67
◆产品类型	Monoclonal
◆基因全称	activating transcription factor 4 (tax-responsive enhancer element B67)
◆经测试应用	WB, IHC
◆宿主	Rabbit
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human ATF4
◆理论分子量	39 kDa
◆实际分子量	45-55 kDa
◆Gene ID (NCBI)	468
◆GenBank Accession Number	BC022088
◆UNIPROT	P18848
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。