

Alix Rabbit Monoclonal antibody

【货号】 RM20259

【规格】 50ul/100ul

【产品信息】

◆产品名称	Alix Rabbit Monoclonal antibody
◆规格	50ul/100ul
◆别名	PDCD6IP, AIP1, ALG-2-interacting protein 1, ALG-2-interacting protein X, DRIP4
◆产品类型	Monoclonal
◆基因全称	programmed cell death 6 interacting protein
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human Alix
◆理论分子量	96 kDa
◆实际分子量	80-100 kDa
◆Gene ID (NCBI)	10015
◆GenBank Accession Number	BC020066
◆UNIPROT	Q8WUM4
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。