

TDP-43(N-terminal) Rabbit Monoclonal antibody

【货号】 RM20271

【规格】 50ul/100ul

【产品信息】

◆产品名称	TDP-43(N-terminal) Rabbit Monoclonal antibody
◆规格	50ul/100ul
◆别名	TDP-43 (for IF/FC), TARDBP, TDP 43, TDP43, TDP-43
◆产品类型	Monoclonal
◆基因全称	TAR DNA binding protein
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human TDP-43
◆理论分子量	43 kDa
◆实际分子量	40 kDa
◆Gene ID (NCBI)	23435
◆GenBank Accession Number	BC001487
◆UNIPROT	Q13148
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:10000-1:50000 for western blotting (WB) 1:100-1:400 for Immunofluorescence (IF/ICC) 1:400-1:1600 for Immunofluorescence (IF-P) 1:800-1:3200 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。