

UQCRC2 Rabbit Monoclonal antibody

【货号】 RM20272

【规格】 50ul/100ul

【产品信息】

◆产品名称	UQCRC2 Rabbit Monoclonal antibody
◆规格	50ul/100ul
◆别名	240454D9, Complex III subunit 2, Core protein II, Cytochrome b-c1 complex subunit 2, mitochondrial, QCR2
◆产品类型	Monoclonal
◆基因全称	ubiquinol-cytochrome c reductase core protein II
◆经测试应用	WB, IHC, IF-P
◆宿主	Rabbit
◆种属反应性	Human, Mouse
◆免疫原	A synthesized peptide derived from human UQCRC2
◆理论分子量	48 kDa
◆实际分子量	48 kDa
◆Gene ID (NCBI)	7385
◆GenBank Accession Number	BC000484
◆UNIPROT	P22695
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF-P) 1:25-1:250 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。