

Anti-alpha tubulin antibody

【货号】 AP10001

【规格】 50ul

【产品信息】

◆产品名称	Anti-alpha tubulin antibody
◆规格	50ul
◆别名	alpha Tubulin, Alpha tubulin ubiquitous, TUBA1B, tubulin, Tubulin alpha, Tubulin alpha 1B chain, Tubulin alpha ubiquitous chain, Tubulin K alpha 1, tubulin, alpha 1b, α -tubulin
◆产品类型	内参抗体
◆基因全称	tubulin, alpha 1b
◆经测试应用	WB, IF/ICC, IHC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human alpha tubulin
◆克隆类型	Polyclonal
◆理论分子量	50 kDa
◆实际分子量	50-55 kDa
◆Gene ID (NCBI)	10376
◆GenBank Accession Number	BC009314
◆UNIPROT	P68363
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:20000-1:100000 for western blotting (WB) 1:500-1:2000 for Immunofluorescence (IF) 1:200-1:1000 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。