

HuR/ELAVL1 Rabbit Polyclonal antibody

【货号】 AP20255

【规格】 50ul/100ul

【产品信息】

◆产品名称	HuR/ELAVL1 Rabbit Polyclonal antibody
◆规格	50ul/100ul
◆别名	ELAVL1, ELAV1, ELAVL 1, ELAV-like protein 1, Hu antigen R
◆产品类型	Polyclonal
◆基因全称	ELAV (embryonic lethal, abnormal vision, Drosophila)-like 1 (Hu antigen R)
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human HuR/ELAVL1
◆理论分子量	36 kDa
◆实际分子量	30-34 kDa
◆Gene ID (NCBI)	1994
◆GenBank Accession Number	BC003376
◆UNIPROT	Q15717
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:25-1:250 for Immunohistochemistry (IHC) 1:25-1:250 for Immunofluorescence (IF/ICC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。