

NCOR2 Rabbit Polyclonal antibody

【货号】 AP20257

【规格】 50ul/100ul

【产品信息】

◆产品名称	NCOR2 Rabbit Polyclonal antibody
◆规格	50ul/100ul
◆别名	CTG repeat protein 26, CTG26, N CoR2, NCOR2, Nuclear receptor corepressor 2, SMAP270, SMRT, SMRTE, SMRTE tau, T3 receptor associating factor, TNRC14, TRAC, TRAC 1, TRAC1
◆产品类型	Polyclonal
◆基因全称	nuclear receptor co-repressor 2
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Rabbit
◆种属反应性	Human
◆免疫原	A synthesized peptide derived from human NCOR2
◆理论分子量	274 kDa
◆实际分子量	230-250 kDa
◆Gene ID (NCBI)	9612
◆GenBank Accession Number	NM_006312
◆UNIPROT	Q9Y618
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:8000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:200-1:800 for Immunohistochemistry (IHC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。